Vertebroplastie Kanüle

zur Zementapplikation

Vertebroplasty Canula

for Bone Cement Application

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INSTRUCTIONS FOR USE

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CE 0123

SOMATEX®

Medical Technologies GmbH
Made in Germany

INSTRUCTIONS FOR USE

ISomatex® - Vertebroplasty Canula serves for percutaneous filling of benign/malignant bone lesions and/or osteoprotic vertebral bodies with bone cement under radiological control. This development corresponds to the very latest available technology, and is characterised by the greatest possible convenience quality and reliability for both the patient and the doctor.

Somatex® - VERTEBROPLASTY-CANULA is not designed for bone marrow biopsies. It is forcefully indicated to respect the processing information in IFU's of used bone cement

Before using pay attention to following steps:

1. Assure that packaging is not opened and/or damaged – check date of shelf life.
2. Checking the needle size.
3. Open the packaging.

Note:

This product has been explicitly designed and manufactured for single use only. The company Somatex® will not be held responsible in any way for the consequences in case of resterilisation

Steps of application:

4. Disinfection of puncture point / covering area by using sterile drape
5. Skin marking of affected vertebral body
6. To facilitate needle insertion make an skin incision by using a scalpel (Fig.11).
7. Plead the needle till resistance of bone tissue becomes perceptible. Following fixing the needle with carefully hammer strokes.
8. Checking of right needle tip position and canula angulations under CT/DSA-Control. Mark at the canula hub shows position of needle tip. This eases orientation in case of necessary corrections of needle direction.
9. After finding right direction drive the canula toward the inside of vertebral body by further careful hammer strokes. During this procedure needle direction has to be controlled under used imaging method.
10. After reaching the target zone and placement of the canula in the right vertebral body segment careful withdraw of canula stilet follows.
11. Preparation of bone cement (Pay attention to instructions for use of bone cement)

12. Connecting syringe with sucked up bone cement to the canula hub and injection of the bone cement into the affected vertebral body.
13. After complete cement injection reinsertion of removed stilet into the canula. This step serves for avoidance of building a cement sting. Please consider that the reinsertion causes a further injection of cement into the vertebral body which remained inside the canula cavity. Corresponding to therapy plan it is mindful that the remained amount of cement is useful for stabilisation of the vertebral body. The below mentioned table shows the volume dependent to used canula size.
14. Up to the point till cement gets hardened (working time see instructions of use for bone cement) turning the canula in minimum for 360° in time steps between 30 and 60 seconds. This avoids a fixing of the canula by the cement inside the vertebral body.
15. After complete hardening of cement withdrawal the canula through pulling and simultaneous turning of the canula
16. Take care for puncture point

Indications:

- Painful bone destructions caused by malignant diseases (osteolytic Metastasis, Plasmocysts)
- Painful bone destructions caused by benign diseases (Osteoporosis, Haemangioma)

Contraindications:

For use of Somatex® - VERTEBROPLASTY-CANULA all contraindications of the particular field of application valid as are known for other needles of this type , given that they are used according to state-of-the-art medical technology. Pay also attention to contraindications which are given in the instructions for use relating individually to each used bone cement

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DECLARATION OF CONFORMITY

In the sense of EEC guideline on Medical Device Directive 93/42/EEC

Manufacturer: **SOMATEX® Medical Technologies GmbH**

Address: Rheinstraße 7d, D-14513 Teltow

Product: **Vertebroplastie-Set**
(REF.: 180770, ..775, ..780, ..785, ..790, ..795, ..800, ..805, ..810)

Classification
(MDD, Annex XI): **IIa**

**We herewith declare that the above mentioned product(s)
corresponds to the following guideline:**

Annex I of EEC-guideline for Medical Devices Directive 93/42/EEC.

Notified Body: TÜV Product Service GmbH, Ridlerstr. 31,
80339 Munich

Certificate: G1 04 08 31459 011

*Date CE mark
was first affixed:* 1st of March 2002

Teltow, 3rd of January 2005

REORDER INFORMATION

| Reorder Number (REF) | Name | Vertebroplasty Canula [Diameter/Length] | Kit components |
|-------------------------|--------------------|--|---|
| 180 795 | Apart-Standard-Kit | 3,0 (11G) / 100mm | 1 Vertebroplasty Canula 1 Polycarbonate Syringe 3ml 1 Power Clip |
| 180 800 | | 3,0 (11G) / 125mm | |
| 180 805 | | 3,0 (11G) / 150mm | |
| 180 810 | | 2,5 (11G) / 125mm | |
| 180 780 | Dual-Kit | 3,0 (11G) / 100mm | 2 Vertebroplasty Canula 2 Polycarbonate Syringe 3ml 2 Power Clip |
| 180 785 | | 3,0 (11G) / 125mm | |
| 180 790 | | 3,0 (11G) / 150mm | |
| 180 770 | Dual-Histo-Kit | 3,0 (11G) / 100mm | 2 Vertebroplasty Canula 1 Biopsy Canula [2,10 (14G) / 150mm] 2 Polycarbonate Syringe 3ml 2 Power Clip |
| 180 775 | | 3,0 (11G) / 150mm | 2 Vertebroplasty Canula 1 Biopsy Canula [2,10 (14G) / 200mm] 2 Polycarbonate Syringe 3ml 2 Power Clip |

*** Every kit contains 1 Plastic Bowl, 1 Spatula, 1 Clamp

SOMATEX®
Medical Technologies GmbH
Made in Germany

Rheinstraße 7d | D-14513 Teltow
 Tel +49 (0) 3328 30 76 - 0 | Fax +49 (0) 3328 30 76 - 99
www.somatex.com | info@somatex.com