OP-Technik Surgical technique

krypton


multisegmental internal fixator

ulrich
medical

OP-Technik
Surgical technique

krypton

Contents

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The small pictures at the bottom of the page show the chronological application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

CE 0123

Art. Nr./Art. No.	Art. Nr./Art. No.	

The spinal fixator krypton is the universal system for stable fixation and reduction of the thoracic and lumbar spine. The constrained assembly at set angles is achieved by means of transpedicular fixation. Mono- and multisegmental instrumentations are possible with the krypton system.

The greatest advantage of the system is the controlled, three-dimensional correction by means of the repositor. The adaptors that connect screw and repositor are used as lever arms. The reduction instrument is outside of the situs and only the reduction arms reach into the wound. This provides a maximum overview without restriction during the reduction, in held position as well as for the exploration for the fusion.

The special feature of the krypton pedicle screw is the connecting cone which is at 45° in regards to the axis of the screw. Thus converging pedicle placement is supported. The axis of the cone faces upwards (top-loading) and makes all further steps easy. The cone is responsible for the very high stability of the screw-rod connection. By this, the screw is fixed to the rod in an optimal way and a high constrained connection at variable angles is achieved. The necessity to use cross links is largely reduced due to the high rotation resistance of the rod-screw interface. Special sacral screws improve the sacral anchorage of the implants. The distance and the angle between the krypton pedicle screws may be freely determined according to the anatomical requirements.

This surgical technique describes the implants and instruments as well as the individual steps for applying krypton. The surgical technique is not sufficient as the sole basis for the successful application of the fixator krypton. It is recommended to study and to learn the operating technique with and from an experienced surgeon. Please note the instructions for reprocessing treatment for implants and instruments.

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Indications

Misalignment, instability, destruction (e.g. fracture, osteolysis), deformity (e.g. kyphosis, spondylolisthesis, scoliosis) and degenerative disorders.

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Contraindications

krypton must not be used in cases of acute or chronic general infection and known intolerance to the material used.
Also krypton must not be used in patients with previous bone diseases such as severe osteoporosis, osteopenia, or osteomyelitis.
Contraindication is a poor medical or psychological general state of the patient.

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General information

The individual anatomy of the patients must allow for the application of pedicle screws. The features and quality of human bone limit the size and resistance of the implant.

Indications and contraindications

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

Spinal Injuries Centre, Stoke Mandeville, UK

Berstungs-Spaltbruch L1, Reposition und Transfixation Th12-L2 mit krypton, interkorporelle Fusion Th12-L1
Burst split fracture L1, reduction and transfixation Th12-L2 with krypton, intersomatic fusion Th12-L1

m., 59 J., pulmonales Karzinom, dorsale Instrumentation Th5-Th7 mit krypton, ventrale Stabilisierung mit ADD, jeweils mit Arthrodesis
m, 59 yrs, pulmonary carcinoma, posterior instrumentation Th5-Th7 with krypton, anterior stabilization with ADD, all with arthrodesis

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krypton pedicle screw

The krypton pedicle screw's head is at an angle of 45° to the axis of the screw. The screw head is equipped with a cone and the thread is self-tapping. The diameter is 5.0 mm, 6.0 mm and 7.5 mm. Screw lengths are available in 5 mm increments: screws Ø 5.0 mm from 35 to 65 mm, screws Ø 6.0 mm and 7.5 mm from 40 to 75 mm.

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krypton sacral screw

The krypton sacral screw is used for an additional anchorage in the sacrum and is positioned divergently in the lateral mass. The diameter is 6.0 mm and the screws are available from 40 to 75 mm length. The axis of the cone is congruent with the axis of the sacral screw.

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Note that the lengths of krypton pedicle screws and sacral screws are determined differently. The picture shows a krypton pedicle and sacral screw of the same length.

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krypton connector

In unconstrained position the connector rotates 360° around the rod as well as around the cone. Therefore a three-dimensional adjustment is possible. Two connectors of different length are available: 19 mm (green) and 21 mm (red).

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krypton locking screw

The locking screw is inserted into the cone. Thus a rotation resistant connection is achieved between connector and rod as well as between connector and screw. The assembly is easy due to the threadless portion of the locking screw.

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krypton rod

The diameter of the krypton rod is 6.25 mm. Rods are available in various lengths from 30 – 500 mm.

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Cross link

- 1) X-plate for crosslink, length 50 and 65 mm
- 2) Hook for crosslink (with a break off notch in the threaded portion)
- 3) Nut for crosslink
- 4) Locking plate for crosslink

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Screw tightener for
krypton pedicle screw

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Drill guide for tissue
protector,
Ø 3.5 mm, 4.5 mm,
6.0 mm

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Probe with spherical
handle and depth
markings, delicate

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Wrench for krypton sacral
screw, hex 9 mm

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Torque wrench with coupling, T-handled (6 and 9 Nm), hex 5 mm

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Tap for krypton pedicle screw, Ø 5.0 mm, 6.0 mm, 7.5 mm

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Repositor incl. 4, 6, 8 or 10 arms.
The repositor consists of the following instruments:
CS 2840-1 to CS 2840-16

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Locking nut for adaptor to
attach to traction device
(CS 2840-9)

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Traction device for
reduction scaffold

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Rod holder for smooth rod,
straight, length 190 mm

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
Awl acc. to Perthes,
T-handle, length 220 mm

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
Exposure of the spine

On the lateral side the posterior segment is exposed including the joints. Using the awl (UT 1850-23) the cortical bone is penetrated under C-arm control. The drilling angle is determined.

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
Preparing the pedicle

The position for the screw insertion must be intraosseous. The bone may be penetrated using a drill bit (CS 2822). Drill through the drill guide (CS 2821-2, 3 or 4) which is attached to the tissue protector (CS 2821-1). Use the following drill bit: drill bit Ø 3.5 mm for 5.0 mm screws, drill bit Ø 4.5 mm for 6.0 mm screws, drill bit Ø 6.0 mm for 7.5 mm screws. Limited by the anatomy, the direction of drilling should be as convergent as possible in order to achieve a high primary stability and to facilitate the implant assembly.

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Preparing the pedicle

As an alternative to the drilling technique the drilling hole may be prepared using an awl or probe with spherical handle (CS 2826). The calibrator on the probe helps to determine the appropriate screw length.

Important:

Prepare the pedicle under C-arm control.

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Exploration of the drill hole

Use the ball tip probe (UL 8522-27) to make sure the pedicle is intact.

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Measuring the drill hole

The screw length may be determined by means of the depth probe (CS 2831).

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

Inserting and securing the adaptors

The krypton screw of the appropriate length is connected to the adaptor (CS 2840-1). The adaptor holds the screw head firmly.


With the screwdriver, 2.5 mm (CS 2813), the pedicle screw is secured in the adaptor.

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

Implanting the pedicle screw

The krypton pedicle screw is implanted by means of the T-wrench for adaptor (CS 2840-2). The cone of the pedicle screw looks inwards. The tap (CS 2832) may be used to prepare the screw implantation, if necessary. The adaptor may also be used as a lever arm and serves as an instrument that holds back the muscles.

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
Assembly of the unconstrained spindle arms

The unconstrained spindle arms (CS 2840-5-1) are connected to the spindle (CS 2840-16, -3 or -4) by means of the adjustor for distraction (CS 2840-5-2). The distance between the arms may be adjusted.

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Assembly of the constrained spindle arm

The constrained spindle arm (CS 2840-5-3) is attached and locked by means of the locker (CS 2840-5-4).

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
Mounting the reduction arms

The reduction arms (CS 2840-6) are connected to the adaptors (CS 2840-1). The muscle retractor (CS 2840-13) serves as a guide so as to facilitate the application of the reduction arms.

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Fixing the reduction arms

The reduction arms are locked with the locking nut (CS 2840-7-1 or -2) by means of a wrench (CS 2840-12). The locking nut CS 2840-7-1 is used if the traction device is not necessary for the reduction.

	CS 2840-6		CS 2840-1		CS 2840-13		CS 2840-7-1/-2		CS 2840-12		CS 2840-16,-3,-4
--	------------------	--	------------------	--	-------------------	--	-----------------------	--	-------------------	--	-------------------------

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Completion of the repositor

The reduction arms are locked to the spindle arms using the locker (CS 2840-5-4). The spindle (CS 2840-16, -3 or -4) with the complete arms is attached to the pedicle screws. The number of the arms may freely be determined. Modular design and length adjustability of the reduction arms make the adaptation to the individual anatomy possible.

Positionierter Reposeur

Einseitige Montage des Reposeurs.

Repositor mounted

Attachment of the repositor on one side.

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
In cases of spondylolisthesis the traction device (CS 2840-9) and the reduction scaffold (CS 2840-8) are used for correction. Usually the hinges of the reduction arms do not need to be locked when an spondylolisthesis is reduced. The locking nut for adaptor to attach to CS 2840-9 has to be used (CS 2840-7-2).

To correct the slippage, the traction device is attached to the vertebral body that has slipped anteriorly. The reduction scaffold is supported by the adjacent segment. Turning the traction device (counter clockwise) reduces the anteriorly slipped vertebral body posteriorly.

The counter piece of the locking nut (CS 2840-7-2) for the tip of the traction device rotates freely and adjusts to the direction of the traction. In order to realign a second vertebral body an additional traction device may be used. Segmental rotation deformities are corrected by means of the traction device and the reduction scaffold as well. Two reduction scaffolds are attached on the posterior side of the involved vertebral bodies and connected to the traction device. Reducing asymmetrically the vertebral bodies derotates the spine.

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In cases of fractures the anatomical profile must be restored and fixed. Lordosation by means of the adaptors that serve as lever arms and simultaneous compression by means of the reduction arms, the deformity is corrected. In held correction the anterior void in the cancellous bone may be filled from the posterior approach if an additional anterior approach is to be avoided.

In held correction the free approach to the segment involved allows for an intralaminar fenestration to decompress the spinal canal or to reconstitute the anterior column.

After correction of fractures the reduction arms have to be fixed with the screwdriver, 3.5 mm (CS 2812).

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Determining the length and contour of the rod

The phantom rod (CS 2819) serves to determine the length and the contour of the rod.

Alternatively the length of the rod may be determined by the measuring caliper (CS 5788) if the rod is to be bent slightly or not at all.

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Contouring the rod

Using the rod bending forceps acc. to Lepine (UT 1639-29), the rod will be brought into the form specified by the phantom rod.

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Assembly of the connectors

If rods are contoured, the connectors can be advanced on the contoured portions of the rods by means of the connector forceps (CS 2818). Two sizes of connectors 19 mm (green) and 21 mm (red) are available.

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
Assembly of the rods

The krypton rod including the pre-assembled connectors may be positioned properly with the rod holder for smooth rod (CS 2781).

The connector will be adjusted to the cone of the pedicle screw. Turning the pedicle screw with the T-wrench (CS 2840-2) adjusts the head of the screw cone to the opening at the connector.

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Applying the connectors

Using the persuader for krypton connector (CS 2825-1 and CS 2825-2) the connectors may be positioned properly so as to match the position of the cone. Introducing the persuader (CS 2825-1) into the cone through the connector. Press the connector onto the cone by means of the sleeve (CS 2825-2) now. Although the persuader is removed the connector remains locked to the cone.

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Fixing the rods


The krypton locking screws are tightened in the cone using the screwdriver (CS 2812). If the situs is very narrow a cardan screwdriver (CS 2814), should be used. The final tightening is performed after all the locking screws have been inserted.

The torque wrench (CS 2828-20, -21) is used to tighten the locking screws to 9 Nm. The finally locking with the torque wrench is done before dismantling the repositor.

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After dismantling the repositor the adaptors are removed using the remover (CS 2830).

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Alternatively the ball tip rod (CS 2840-15) could be used to remove the adapters.

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

The sacral purchase can be augmented in case of poor bone quality. The pedicle screws are equipped with rods that exceed approximately 30 mm caudally. The rods are contoured, pointing medially. The exceeding portion is equipped with an additional connector.

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
Drilling for sacral anchorage

Using the tissue protector (CS 2821-1) and the drill guide (CS 2821-3) drill through the attached connector. The probe with spherical handle (CS 2826) may be used as an alternative to the drill (CS 2822-2). Determine the length of the sacral screw using the depth probe (CS 2831) through the connector.

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
Implanting the sacral screw

The connector is rotated back in order to screw in the sacral screw divergently in the lateral mass, using the screwdriver, 3.5 mm (CS 2812).

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Fixing the rods

The krypton locking screw is introduced into the sacral rod through the connector, using the screwdriver, 3.5 mm (CS 2812). The sacral screw is secured against rotation by means of the wrench for sacral screw (CS 2827).

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Complete assembly with additional sacral screws

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

[A] Insert the hooks for cross link (CS 2809-3). Hooks can be placed laterally (picture) or medially.

[B] Attach the X-plate for cross link (CS 2809-1/2).

[C] Attach the locking plate (CS 2809-5) and fix the nut for cross link (CS 2809-4) with the socket wrench for nut (CS 2815).

[D] Trim the rod of the hooks with the socket wrench for nut.

[E] The torque wrench (CS 2828-20, CS 2828-22) is used to tighten the nut to 6 Nm.

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Complete assembly with cross link stabilization

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Removal of connector

To disconnect the cone, the remover for krypton connector (CS 2824-2) is used.

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Removal of pedicle screws

Removal of the pedicle screws is possible with the screw tightener (CS 2811).

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
Layer 3 and tray for krypton instruments

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Implants

Product number

krypton pedicle screw, Ø 5.0 mm, length 35 mm	CS 2800-5-35T
krypton pedicle screw, Ø 5.0 mm, length 40 mm	CS 2800-5-40T
krypton pedicle screw, Ø 5.0 mm, length 45 mm	CS 2800-5-45T
krypton pedicle screw, Ø 5.0 mm, length 50 mm	CS 2800-5-50T
krypton pedicle screw, Ø 5.0 mm, length 55 mm	CS 2800-5-55T
krypton pedicle screw, Ø 5.0 mm, length 60 mm	CS 2800-5-60T
krypton pedicle screw, Ø 5.0 mm, length 65 mm	CS 2800-5-65T
krypton pedicle screw, Ø 6.0 mm, length 35 mm	CS 2800-6-35T
krypton pedicle screw, Ø 6.0 mm, length 40 mm	CS 2800-6-40T
krypton pedicle screw, Ø 6.0 mm, length 45 mm	CS 2800-6-45T
krypton pedicle screw, Ø 6.0 mm, length 50 mm	CS 2800-6-50T
krypton pedicle screw, Ø 6.0 mm, length 55 mm	CS 2800-6-55T
krypton pedicle screw, Ø 6.0 mm, length 60 mm	CS 2800-6-60T
krypton pedicle screw, Ø 6.0 mm, length 65 mm	CS 2800-6-65T
krypton pedicle screw, Ø 6.0 mm, length 70 mm	CS 2800-6-70T
krypton pedicle screw, Ø 6.0 mm, length 75 mm	CS 2800-6-75T
krypton pedicle screw, Ø 7.5 mm, length 40 mm	CS 2800-7-40T
krypton pedicle screw, Ø 7.5 mm, length 45 mm	CS 2800-7-45T
krypton pedicle screw, Ø 7.5 mm, length 50 mm	CS 2800-7-50T
krypton pedicle screw, Ø 7.5 mm, length 55 mm	CS 2800-7-55T
krypton pedicle screw, Ø 7.5 mm, length 60 mm	CS 2800-7-60T
krypton pedicle screw, Ø 7.5 mm, length 65 mm	CS 2800-7-65T
krypton pedicle screw, Ø 7.5 mm, length 70 mm	CS 2800-7-70T
krypton pedicle screw, Ø 7.5 mm, length 75 mm	CS 2800-7-75T
krypton sacral screw, Ø 6.0 mm, length 40 mm	CS 2801-6-40T
krypton sacral screw, Ø 6.0 mm, length 45 mm	CS 2801-6-45T
krypton sacral screw, Ø 6.0 mm, length 50 mm	CS 2801-6-50T
krypton sacral screw, Ø 6.0 mm, length 55 mm	CS 2801-6-55T
krypton sacral screw, Ø 6.0 mm, length 60 mm	CS 2801-6-60T
krypton sacral screw, Ø 6.0 mm, length 65 mm	CS 2801-6-65T
krypton sacral screw, Ø 6.0 mm, length 70 mm	CS 2801-6-70T
krypton sacral screw, Ø 6.0 mm, length 75 mm	CS 2801-6-75T
krypton sacral screw, Ø 6.0 mm, length 85 mm	CS 2801-6-85T
krypton sacral screw, Ø 6.0 mm, length 95 mm	CS 2801-6-95T
krypton connector, green, length 19 mm	CS 2805-19T
krypton connector, red, length 21 mm	CS 2805-21T
krypton locking screw	CS 2806T
krypton rod, Ø 6.25 mm, length 30 mm	CS 2808-30T
krypton rod, Ø 6.25 mm, length 35 mm	CS 2808-35T
krypton rod, Ø 6.25 mm, length 40 mm	CS 2808-40T
krypton rod, Ø 6.25 mm, length 45 mm	CS 2808-45T
krypton rod, Ø 6.25 mm, length 50 mm	CS 2808-50T
krypton rod, Ø 6.25 mm, length 55 mm	CS 2808-55T
krypton rod, Ø 6.25 mm, length 60 mm	CS 2808-60T
krypton rod, Ø 6.25 mm, length 65 mm	CS 2808-65T
krypton rod, Ø 6.25 mm, length 70 mm	CS 2808-70T
krypton rod, Ø 6.25 mm, length 75 mm	CS 2808-75T
krypton rod, Ø 6.25 mm, length 80 mm	CS 2808-80T
krypton rod, Ø 6.25 mm, length 85 mm	CS 2808-85T
krypton rod, Ø 6.25 mm, length 90 mm	CS 2808-90T
krypton rod, Ø 6.25 mm, length 95 mm	CS 2808-95T
krypton rod, Ø 6.25 mm, length 100 mm	CS 2808-100T
krypton rod, Ø 6.25 mm, length 108 mm	CS 2808-108T
krypton rod, Ø 6.25 mm, length 120 mm	CS 2808-120T
krypton rod, Ø 6.25 mm, length 132 mm	CS 2808-132T
krypton rod, Ø 6.25 mm, length 144 mm	CS 2808-144T
krypton rod, Ø 6.25 mm, length 156 mm	CS 2808-156T
krypton rod, Ø 6.25 mm, length 168 mm	CS 2808-168T
krypton rod, Ø 6.25 mm, length 180 mm	CS 2808-180T
krypton rod, Ø 6.25 mm, length 192 mm	CS 2808-192T
krypton rod, Ø 6.25 mm, length 204 mm	CS 2808-204T
krypton rod, Ø 6.25 mm, length 216 mm	CS 2808-216T
krypton rod, Ø 6.25 mm, length 375 mm	CS 2808-375T
krypton rod, Ø 6.25 mm, length 500 mm	CS 2808-500T

Implants

	Product number
X-plate for crosslink, length 50 mm	CS 2809-1T
X-plate for crosslink, length 65 mm	CS 2809-2T
Hook for crosslink	CS 2809-3T
Nut for crosslink	CS 2809-4T
Locking plate for crosslink	CS 2809-5T

Instruments

	Product number
Screw tightener for krypton pedicle screw	CS 2811
Screwdriver, hex 3.5 mm	CS 2812
Screwdriver, hex 2.5 mm for CS 2840-1	CS 2813
Cardan screwdriver	CS 2814
Socket wrench for nut (CS 2809-4T), hex 8 mm	CS 2815
Connector forceps	CS 2818
Phantom rod, length 150 mm	CS 2819-150
Phantom rod, length 400 mm	CS 2819-400
Tissue protector	CS 2821-1
Drill guide for tissue protector, Ø 3.5 mm	CS 2821-2
Drill guide for tissue protector, Ø 4.5 mm	CS 2821-3
Drill guide for tissue protector, Ø 6.0 mm	CS 2821-4
Drill bit, Ø 3.5 mm	CS 2822-1
Drill bit, Ø 4.5 mm	CS 2822-2
Drill bit, Ø 6.0 mm	CS 2822-3
AO drill bit, Ø 3.5 mm	CS 2822-4
AO drill bit, Ø 4.5 mm	CS 2822-5
AO drill bit, Ø 6.0 mm	CS 2822-6
Remover for krypton connector	CS 2824-2
Persuador for krypton connector	CS 2825-1
Sleeve for CS 2825-1	CS 2825-2
Probe with spherical handle and depth markings, delicate	CS 2826
Wrench for krypton sacral screw, hex 9 mm	CS 2827
Torque wrench with coupling, T-handled (6 and 9 Nm), hex 5 mm	CS 2828-20
Hex bit, 3.5 mm for CS 2828-20	CS 2828-21
Hex socket, 8.0 mm for CS 2828-20	CS 2828-22
Remover for adaptor	CS 2830
Depth probe	CS 2831
Tap for krypton pedicle screw, Ø 5.0 mm	CS 2832-1
Tap for krypton pedicle screw, Ø 6.0 mm	CS 2832-2
Tap for krypton pedicle screw, Ø 7.5 mm	CS 2832-3
Adaptor	CS 2840-1
T-wrench for adaptor	CS 2840-2
Spindle, length 500 mm	CS 2840-3
Spindle, length 300 mm	CS 2840-4
Spindle arm, unconstrained	CS 2840-5-1
Adjustor for distraction	CS 2840-5-2
Spindle arm, constrained	CS 2840-5-3
Locker	CS 2840-5-4
Reduction arm	CS 2840-6
Locking nut for adaptor	CS 2840-7-1
Locking nut for adaptor to attach to traction device (CS 2840-9)	CS 2840-7-2
Reduction scaffold	CS 2840-8
Traction device for reduction scaffold	CS 2840-9
Wrench for locking nut	CS 2840-12
Muscle retractor	CS 2840-13
Distraction helper for CS 2840-5-2	CS 2840-14
Ball tip rod, Ø 3.0 mm	CS 2840-15
Spindle, length 170 mm	CS 2840-16
Rod holder for smooth rod, straight, length 190 mm	CS 2781
Measuring caliper	CS 5788
Ball tip probe	UL 8522-27
Rod bending forceps acc. to Lepine, length 290 mm	UT 1639-29
Awl acc. to Perthes, T-handle, length 220 mm	UT 1850-23

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ulrich GmbH & Co. KG