

SURGICAL TECHNIQUE
FIBULAR DISTAL PLATE ANGULARLY STABLE



Fibular Distal Plates

Surgical technique	1–5
Implants	6–7
Instruments	8

Description of the medical device

The implant system consists of a plate and the necessary number of angularly stable screws, or cortical screws.

Fibular distal plate

The plates are designed to be non-side-specific. They are anatomically shaped according to the distal end of the fibula.

The plates come in four sizes: 4–10-hole, increasing by 2 holes, or 4 lengths, 57–105 mm increasing by 16 mm.

The holes in the plate are perpendicular to the axis of the fibular.

The plates come in titanium and stainless steel versions.



Angularly stable screws

The distal fibular plate is fixed using Ø 3.5 mm and Ø 3.5/2.7 mm angularly stable screws.

Ø 3.5 mm angularly stable screws come in lengths of 8–50 mm increasing by 2 mm and 50–75 mm by 5 mm.

Ø 3.5/2.7 mm angularly stable screws come in lengths of 10–50 mm increasing by 2 mm.




Cortical screws

Cortical HA 3.5 mm screws are intended for the fibular plate.

Indications

The fibular distal plate is intended for angularly stable osteosynthesis of the distal end of the fibular bone, primarily for treating fracture types 44 – A1, A2, B1, B3, C1, C2, C3 according to Müller.



This brochure only serves as an illustrative guide for fibular distal plates and the instrumentation. The aim of the brochure is to give physicians and suture nurses a quick guide and the correct composition in the use of the instruments and implant in order to achieve the best operational outcome. If you have any questions, please contact your MEDIN, a.s. salesperson.

Surgical technique

1. Operational planning

Before the osteosynthesis we recommend that you go through the surgical procedure. Pre-select the appropriate implants. Compare the X-Ray photos of both limbs and then specify the repositioning of the fragments back into the anatomical position. With regards to the soft tissues, we recommend planning the location of the supporting implants (K-wires, spongy screws) used during the transfixation of fragments and the main implants.

2. Patient position

Place the patient on the operating table on their back so that the limb to be operated upon is in neutral rotation. This provides support (fig. 2). We recommend placing the patient so that they can be checked using two X-Ray projections perpendicular to one another.

3. Operational approach

Make an incision laterally along the fibular bone. The incision length is determined by the length of the plate. If appropriate, the plate can be introduced using mini-invasive techniques, a short incision around 50 mm. During mini-invasive introduction, use a raspator to create a space for introducing the plates. Only use mini-invasive techniques to fix simple fractures, which do not require a visual check of the fragments' repositioning. Keep in mind that there may be arteries, veins and nerves at the incision site. Damage to these tissues, nerves especially, can have lasting consequences for the patient.

4. Repositioning

Reposition the fragments into the original anatomical position. The repositioned fragments are then fixed with at least two K-wires, so as not to obstruct the introduction of the plate. Check the repositioning by X-Ray. Keep in mind that poor repositioning may adversely affect the fragments' correct working and adhesion.

5. Introducing the plate

Introduce a lockable aiming sleeve to the plate, this improves handling of the plate during introduction (fig. 3). Introduce the plate and fix it to the oval aperture using a 3.5 HA cortical screw. Drill a hole Ø 2.9 mm (fig. 4) and introduce the cortical screw with a 2.5 mm hexagon screwdriver (fig. 5). This enables the plate's position to be corrected in the direction of the fibula axis.




Fig. 2



Fig. 3



Fig. 4



Fig. 5

To fix the plate use a K-wire \varnothing 1.5 mm in the chosen hole through a lockable sleeve and an aiming sleeve (fig. 6).










Fig. 6

AN OVERVIEW OF THE TOOLS:

	Lockable aiming sleeve 2.9 mm; 60 mm	marked green
	Aiming sleeve for 1.5 mm wire; 75 mm	marked green/pink
	Guide rod 1.5 x 300 mm	marked in sieve – pink
	Aiming sleeve for 2.0 mm wires; 75 mm	marked green/purple
	K-wire with shank; 2.0 x 300 mm	marked in sieve – purple
	Drill 2.9 x 190 mm	marked green

6. Introducing a screw

Angularly stable screw \varnothing 3.5 mm

Introduce the lockable sleeve to the selected holes and drill a \varnothing 2.9 mm hole to introduce the angularly stable screw \varnothing 3.5 mm (fig. 7).



Fig. 7

Angularly stable screw \varnothing 3.5/2.7 mm

To introduce a stable \varnothing 3.5/2.7 mm screw introduce a \varnothing 2.9/2 mm aiming sleeve into the lockable sleeve and form a hole using a \varnothing 3.5 mm K-wire (fig. 8).



Fig. 8




Fig. 9



Fig. 10

Remove the sleeve and measure the length of the screw with the depth gauge (Fig. 9). Introduce the screw using the torque limiter, which „clicks“ to indicate the optimum tightening of the screw (Fig. 10).

Implant the required number of screws using the same procedure (Fig. 11).




Fig. 11

7) Introducing the suprasyndesmal screw

If the fibula is unstable in the tibial incisura, it may be indicated by the introduction of a suprasyndesmal screw. In this case, do not primarily fix the plate through the oval aperture, use the sleeves and K-wires. After the osteosynthesis of the fibula, place the fibula in the correct anatomical position in the incisura and temporarily fix with K-wires just above and in parallel with the joint surface. In the proximal part of the oval hole, drill a hole for the HA 3.5 cortical screw – about 2–4 cm above the joint surface of the ankle (fig. 12).



Fig. 12

Measure the required length of the screw and introduce it (fig. 13).



Fig. 13

Introduce the screw tricortically, or, depending on the stability and quality of the bone tissue, quadricortically (fig. 14).

Note: The suprasyndesmal screw is not introduced for tightening, but as a spacer.





Fig. 14

8. Completing the operation

After rinsing out the wound, gradually close the incision. Introduce a suction drain to where the plate has been introduced. The wound is covered with a soft bandage and X-Ray documentation taken.



9. Concluding remarks

- When applied to one patient there should never be a combination of different materials.
- To guarantee the safe use of the implant, MEDIN requires only implants of this company be used. There must not be a combination of implants from other firms.
- The patient must be warned that the implant does not transfer the patient's entire weight. The patient must use a means of support when walking and gradually increase the burden depending on how the muscles are forming in the site of the fracture.
- The implants are intended for single use, for one patient and for one stabilization of the damaged bone. Reuse is forbidden. This fact is mentioned in the package leaflet and concerns all implants.




10. Recommended procedure for extracting the implant

Loosen all the screws and then remove them, this prevents the last screw from rotating along with the splint. In the event of complications, use the instrumentarium to remove the screws.


FIBULAR DISTAL PLATE

FIBULAR DISTAL PLATES ANGULARLY STABLE



SSt	Ti	holes	A	silă
129 78 5180	129 78 5183	4	57 mm	2,2 mm
129 78 5190	129 78 5193	6	73 mm	2,2 mm
129 78 5200	129 78 5203	8	89 mm	2,2 mm
129 78 5210	129 78 5213	10	105 mm	2,2 mm

LOCKING BONE SCREWS 3,5/2,7 mm




129 77 7531
self-tapping
30 mm

SSt	Ti	A
129 77 7431	129 77 7434	10 mm
129 77 7441	129 77 7444	12 mm
129 77 7451	129 77 7454	14 mm
129 77 7461	129 77 7464	16 mm
129 77 7471	129 77 7474	18 mm
129 77 7481	129 77 7484	20 mm
129 77 7491	129 77 7494	22 mm
129 77 7501	129 77 7504	24 mm
129 77 7511	129 77 7514	26 mm
129 77 7521	129 77 7524	28 mm
129 77 7531	129 77 7534	30 mm
129 77 7541	129 77 7544	32 mm
129 77 7551	129 77 7554	34 mm
129 77 7561	129 77 7564	36 mm
129 77 7571	129 77 7574	38 mm
129 77 7921	129 77 7924	40 mm
129 77 7931	129 77 7934	42 mm
129 77 7941	129 77 7944	44 mm
129 77 7951	129 77 7954	46 mm
129 77 7961	129 77 7964	48 mm
129 77 7971	129 77 7974	50 mm

NOTES:
drilled with Ø 2.0 mm wire
introduction by a screwdriver with 2.5 mm socket

NOTES: SSt – stainless steel in accordance with ISO 5832-1 Ti – titanium version, material: Ti6Al4V ELI in accordance with ISO 5832-3

LOCKING BONE SCREWS 3,5 mm




NOTES:
drilled with Ø 2.9 mm drill
introduction by a screwdriver with 2.5 mm socket

SSt	Ti	A
129 77 7021	129 77 7024	8 mm
129 77 7031	129 77 7034	10 mm
129 77 7041	129 77 7044	12 mm
129 77 7051	129 77 7054	14 mm
129 77 7061	129 77 7064	16 mm
129 77 7071	129 77 7074	18 mm
129 77 7081	129 77 7084	20 mm
129 77 7091	129 77 7094	22 mm
129 77 7101	129 77 7104	24 mm
129 77 7111	129 77 7114	26 mm
129 77 7121	129 77 7124	28 mm
129 77 7131	129 77 7134	30 mm
129 77 7141	129 77 7144	32 mm
129 77 7151	129 77 7154	34 mm
129 77 7161	129 77 7164	36 mm
129 77 7171	129 77 7174	38 mm
129 77 7181	129 77 7184	40 mm
129 77 7191	129 77 7194	42 mm
129 77 7201	129 77 7204	44 mm
129 77 7211	129 77 7214	46 mm
129 77 7221	129 77 7224	48 mm
129 77 7231	129 77 7234	50 mm
129 77 7241	129 77 7244	55 mm
129 77 7251	129 77 7254	60 mm
129 77 7261	129 77 7264	65 mm
129 77 7271	129 77 7274	70 mm
129 77 7281	129 77 7284	75 mm



SELF-TAPPING CORTICAL BONE SCREWS – HA 3,5




NOTES:
drilled with Ø 2.9 mm drill
introduction by a screwdriver with 2.5 mm socket


SSt	Ti	A
129 79 5241	129 79 5244	16 mm
129 79 5251	129 79 5254	18 mm
129 79 5261	129 79 5264	20 mm
129 79 5271	129 79 5274	22 mm
129 79 5281	129 79 5284	24 mm
129 79 5291	129 79 5294	26 mm
129 79 5301	129 79 5304	28 mm
129 79 5311	129 79 5314	30 mm
129 79 5321	129 79 5324	32 mm
129 79 5331	129 79 5334	34 mm
129 79 5341	129 79 5344	36 mm
129 79 5351	129 79 5354	38 mm
129 79 5361	129 79 5364	40 mm
129 79 5371	129 79 5374	42 mm
129 79 5441	129 79 5444	44 mm
129 79 5451	129 79 5454	46 mm
129 79 5461	129 79 5464	48 mm
129 79 5391	129 79 5394	50 mm
129 79 5401	129 79 5404	55 mm
129 79 5411	129 79 5414	60 mm
129 79 5421	129 79 5424	65 mm
129 79 5431	129 79 5434	70 mm

FIBULAR DISTAL PLATE


INSTRUMENTS FOR ANGULARLY STABLE PLATES WITH SCREWS 3,5 mm



1	129 09 2571	K-wire with shank; 2,0 mm; 300 mm	pcs
2	129 09 3000	Guide rod 1,5 mm; 300 mm	2
3	129 69 3360	Aiming sleeve lockable 2,9 mm; 60 mm	2
4	129 69 3370	Aiming sleeve for wires 1,5 mm; 75 mm	3
5	129 69 4360	Aiming sleeve for wires 2,0 mm; 75 mm	2
6	129 69 4780	Depth gauge	2
7	129 69 5131	Screwdriver handle	1
8	129 69 5126	Torque limiter 1,5 Nm	1
9	129 69 5231	Screwdriver; hexagon 2,5 mm; 160 mm	1
10	129 79 9981	Drill 2,9 mm; 190 mm	1



SET OF INSTRUMENTS FOR ASP 3,5
139 09 0255
540 × 240 × 50 mm
with instruments



SIEVE FOR INSTRUMENTS FOR ASP 3,5
129 69 4390
540 × 240 × 50 mm
without instruments

© 2011 MEDIN, a.s.; All rights reserved.
This document should be used for commercial purposes of MEDIN, a.s.; the data mentioned in the document has informative character. No part of this document can be copied or published in any form without approval of MEDIN, a.s. The product design may differ from those depicted in these illustrations at the date of issue. Adjustments, made from the reason of further developments of technical parameters, are reserved. Printing and typographical errors are reserved.

REFD